



BIRC5 Antibody, HRP conjugated

Product Code	CSB-PA07249B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15392
Immunogen	Recombinant Human Baculoviral IAP repeat-containing protein 5 protein (1-130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Baculoviral IAP repeat-containing protein 5 (Apoptosis inhibitor 4) (Apoptosis inhibitor survivin), BIRC5, API4 IAP4
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BIRC5