



IL32 Antibody, Biotin conjugated

Product Code	CSB-PA08309D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24001
Immunogen	Recombinant Human Interleukin-32 protein (31-234AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-32 (IL-32) (Natural killer cells protein 4) (Tumor necrosis factor alpha-inducing factor), IL32, NK4 TAIF
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL32
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.