





ANGPTL1 Antibody

Product Code	CSB-PA09689A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95841
Immunogen	Recombinant Human Angiopoietin-related protein 1 protein (24-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiopoietin-related protein 1 (Angiopoietin-3) (ANG-3) (Angiopoietin-like protein 1), ANGPTL1, ANG3 ANGPT3 ARP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ANGPTL1

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA09689A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA09689A0Rb at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.