



CD34 Antibody, Biotin conjugated

Product Code	CSB-PA10509D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28906
Immunogen	Recombinant Human Hematopoietic progenitor cell antigen CD34 protein (124-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hematopoietic progenitor cell antigen CD34 (CD antigen CD34), CD34
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CD34
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.