



CD34 Antibody, Biotin conjugated

Product Code	CSB-PA10509D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28906
Immunogen	Recombinant Human Hematopoietic progenitor cell antigen CD34 protein (124-375AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hematopoietic progenitor cell antigen CD34 (CD antigen CD34), CD34
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CD34