



FCGR2B Antibody, Biotin conjugated

Product Code	CSB-PA10839D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31994
Immunogen	Recombinant Human Low affinity immunoglobulin gamma Fc region receptor II-b protein (12-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor II-b (IgG Fc receptor II-b) (CDw32) (Fc-gamma RII-b) (Fc-gamma-RIIb) (FcRII-b) (CD antigen CD32), FCGR2B, CD32 FCG2 IGFR2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCGR2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.