



A2M Antibody, Biotin conjugated

Product Code	CSB-PA14579D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01023
Immunogen	Recombinant Human Alpha-2-macroglobulin protein (1216-1467AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-2-macroglobulin (Alpha-2-M) (C3 and PZP-like alpha-2-macroglobulin domain-containing protein 5), A2M, CPAMD5
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	A2M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.