





RYR2 Antibody

Product Code	CSB-PA15139A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92736
Immunogen	Recombinant Human Ryanodine receptor 2 protein (901-1192AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ryanodine receptor 2 (RYR-2) (RyR2) (hRYR-2) (Cardiac muscle ryanodine receptor) (Cardiac muscle ryanodine receptor-calcium release channel) (Type 2 ryanodine receptor), RYR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RYR2


Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA15139A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA15139A0Rb at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA15139A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.