



RYR2 Antibody, Biotin conjugated

Product Code	CSB-PA15139D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92736
Immunogen	Recombinant Human Ryanodine receptor 2 protein (901-1192AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ryanodine receptor 2 (RYR-2) (RyR2) (hRYR-2) (Cardiac muscle ryanodine receptor) (Cardiac muscle ryanodine receptor-calcium release channel) (Type 2 ryanodine receptor), RYR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RYR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.