



MLH1 Antibody, Biotin conjugated

Product Code	CSB-PA17009D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40692
Immunogen	Recombinant Human DNA mismatch repair protein Mlh1 protein (350-519AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA mismatch repair protein Mlh1 (MutL protein homolog 1), MLH1, COCA2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MLH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.