



NFATC1 Antibody, Biotin conjugated

Product Code	CSB-PA17827D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95644
Immunogen	Recombinant Human Nuclear factor of activated T-cells, cytoplasmic 1 protein (48-406AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor of activated T-cells, cytoplasmic 1 (NF-ATc1) (NFATc1) (NFAT transcription complex cytosolic component) (NF-ATc) (NFATc), NFATC1, NFAT2 NFATC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NFATC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.