



atpB Antibody, HRP conjugated

Product Code	CSB-PA20469B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19366
Immunogen	Recombinant Arabidopsis thaliana ATP synthase subunit beta, chloroplastic protein (92-393AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP synthase subunit beta, chloroplastic (EC 3.6.3.14) (ATP synthase F1 sector subunit beta) (F-ATPase subunit beta), atpB
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	atpB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.