



PSBO Antibody

Product Code	CSB-PA318952HA01FKI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12359
Immunogen	Recombinant Spinacia oleracea Oxygen-evolving enhancer protein 1, chloroplastic protein (63-232AA)
Raised In	Rabbit
Species Reactivity	Spinacia oleracea
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxygen-evolving enhancer protein 1, chloroplastic (OEE1) (33 kDa subunit of oxygen evolving system of photosystem II) (33 kDa thylakoid membrane protein) (OEC 33 kDa subunit), PSBO
Immunogen Species	Spinacia oleracea
Research Area	Others
Target Names	PSBO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.