



# PSBO Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA318952HC01FKI
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P12359
<b>Immunogen</b>	Recombinant Spinacia oleracea Oxygen-evolving enhancer protein 1, chloroplastic protein (63-232AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Spinacia oleracea
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Oxygen-evolving enhancer protein 1, chloroplastic (OEE1) (33 kDa subunit of oxygen evolving system of photosystem II) (33 kDa thylakoid membrane protein) (OEC 33 kDa subunit), PSBO
<b>Immunogen Species</b>	Spinacia oleracea
<b>Research Area</b>	Others
<b>Target Names</b>	PSBO
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.