



VACWR062 Antibody, FITC conjugated

Product Code	CSB-PA325116LC01VAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21607
Immunogen	Recombinant Vaccinia virus Protein E6 protein (180-429AA)
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E6
Immunogen Species	Vaccinia virus
Research Area	Others
Target Names	VACWR062
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.