



VACWR062 Antibody, Biotin conjugated

Product Code	CSB-PA325116LD01VAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21607
Immunogen	Recombinant Vaccinia virus Protein E6 protein (180-429AA)
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E6
Immunogen Species	Vaccinia virus
Research Area	Others
Target Names	VACWR062