



PRDX2 Antibody, HRP conjugated

Product Code	CSB-PA339235LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32119
Immunogen	Recombinant Human Peroxiredoxin-2 protein (2-65AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxiredoxin-2 (EC 1.11.1.15) (Natural killer cell-enhancing factor B) (NKEF-B) (PRP) (Thiol-specific antioxidant protein) (TSA) (Thioredoxin peroxidase 1) (Thioredoxin-dependent peroxide reductase 1), PRDX2, NKEFB TDPX1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PRDX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.