



NRXN1 Antibody, HRP conjugated

Product Code	CSB-PA347832LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P58400
Immunogen	Recombinant Human Neurexin-1-beta protein (184-363AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurexin-1-beta (Neurexin I-beta), NRXN1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NRXN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.