



Mup6 Antibody, Biotin conjugated

Product Code	CSB-PA360881HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02762
Immunogen	Recombinant Mouse Major urinary protein 6 protein (19-180AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major urinary protein 6 (MUP 6) (Alpha-2U-globulin) (Group 1, BS6) (allergen Mus m 1), Mup6
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Mup6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.