



ISPD Antibody, Biotin conjugated

Product Code	CSB-PA388827LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4D126
Immunogen	Recombinant Human Isoprenoid synthase domain-containing protein (1-264AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	D-ribitol-5-phosphate cytidyltransferase (EC 2.7.7.40) (2-C-methyl-D-erythritol 4-phosphate cytidyltransferase-like protein) (Isoprenoid synthase domain-containing protein) (hISPD), ISPD
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ISPD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.