



MCIDAS Antibody, Biotin conjugated

Product Code	CSB-PA518689LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D6RGH6
Immunogen	Recombinant Human Multicilin protein (3-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multicilin (Multiciliate differentiation and DNA synthesis-associated cell cycle protein) (Mcidas protein) (Protein Idas), MCIDAS, IDAS MCI MCIN
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MCIDAS