



MCIDAS Antibody, Biotin conjugated

Product Code	CSB-PA518689LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D6RGH6
Immunogen	Recombinant Human Multicilin protein (3-110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multicilin (Multiciliate differentiation and DNA synthesis-associated cell cycle protein) (Mcidas protein) (Protein Idas), MCIDAS, IDAS MCI MCIN
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MCIDAS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.