



ADCY1 Antibody, Biotin conjugated

Product Code	CSB-PA600142LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08828
Immunogen	Recombinant Human Adenylate cyclase type 1 protein (416-610AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Adenylate cyclase type 1 (EC 4.6.1.1) (ATP pyrophosphate-lyase 1) (Adenylate cyclase type I) (Adenylyl cyclase 1) (Ca(2+)/calmodulin-activated adenylyl cyclase), ADCY1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADCY1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.