



PTPRS Antibody, HRP conjugated

Product Code	CSB-PA613393LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13332
Immunogen	Recombinant Human Receptor-type tyrosine-protein phosphatase S protein (25-128AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-type tyrosine-protein phosphatase S (R-PTP-S) (EC 3.1.3.48) (Receptor-type tyrosine-protein phosphatase sigma) (R-PTP-sigma), PTPRS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PTPRS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.