





ZRSR2P1 Antibody

Product Code	CSB-PA614898LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15695
Immunogen	Recombinant Human U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 protein (197-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U2 small nuclear ribonucleoprotein auxiliary factor 35 kDa subunit-related protein 1 (CCCH type zinc finger, RNA-binding motif and serine/arginine rich protein 1) (U2(RNU2) small nuclear RNA auxiliary factor 1-like 1), ZRSR1, U2AF1-RS1 U2AF1L1 U2AF1P U2AF1RS1 U2AFBPL
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZRSR2P1

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA614898LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA614898LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.