





PDK4 Antibody

Product Code	CSB-PA614983LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16654
Immunogen	Recombinant Human [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial protein (143-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 4), PDK4, PDHK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDK4

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA614983LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA614983LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.