



PDK4 Antibody, Biotin conjugated

Product Code	CSB-PA614983LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16654
Immunogen	Recombinant Human [Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial protein (143-333AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	[Pyruvate dehydrogenase (acetyl-transferring)] kinase isozyme 4, mitochondrial (EC 2.7.11.2) (Pyruvate dehydrogenase kinase isoform 4), PDK4, PDHK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDK4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.