



DNAJC3 Antibody, FITC conjugated

Product Code	CSB-PA615670LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13217
Immunogen	Recombinant Human DnaJ homolog subfamily C member 3 protein (277-496AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily C member 3 (Endoplasmic reticulum DNA J domain-containing protein 6) (ER-resident protein ERdj6) (ERdj6) (Interferon-induced, double-stranded RNA-activated protein kinase inhibitor) (Protein kinase inhibitor of 58 kDa) (Protein kinase inhibitor p58), DNAJC3, P58IPK PRKRI
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	DNAJC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.