



SQSTM1 Antibody, Biotin conjugated

Product Code	CSB-PA615696LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13501
Immunogen	Recombinant Human Sequestosome-1 protein (16-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sequestosome-1 (EBI3-associated protein of 60 kDa) (EBIAP) (p60) (Phosphotyrosine-independent ligand for the Lck SH2 domain of 62 kDa) (Ubiquitin-binding protein p62), SQSTM1, ORCA OSIL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SQSTM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.