



STIL Antibody, Biotin conjugated

Product Code	CSB-PA618008LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15468
Immunogen	Recombinant Human SCL-interrupting locus protein (388-637AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SCL-interrupting locus protein (TAL-1-interrupting locus protein), STIL, SIL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	STIL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.