



LPAL2 Antibody, Biotin conjugated

Product Code	CSB-PA618087LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16609
Immunogen	Recombinant Human Putative apolipoprotein(a)-like protein 2 protein (52-124AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative apolipoprotein(a)-like protein 2 (Apo(a)-like protein 2) (Lp(a)-liker protein 2) (Apolipoprotein a-related gene C protein) (Apo(a)rg-C), LPAL2, APOARGC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LPAL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.