




ARHGEF5 Antibody

Product Code	CSB-PA618628LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12774
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 5 protein (433-507AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho guanine nucleotide exchange factor 5 (Ephexin-3) (Guanine nucleotide regulatory protein TIM) (Oncogene TIM) (Transforming immortalized mammary oncogene) (p60 TIM), ARHGEF5, TIM
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGEF5

Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA618628LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.