




HIST3H3 Antibody

Product Code	CSB-PA623092HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16695
Immunogen	Recombinant Human Histone H3.1t protein (2-30AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H3.1t (H3/t) (H3t) (H3/g), HIST3H3, H3FT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIST3H3

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA623092HA01HU at dilution of 1: 100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.