



BRD3 Antibody, HRP conjugated

Product Code	CSB-PA624013HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15059
Immunogen	Recombinant Human Bromodomain-containing protein 3 protein (356-506AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bromodomain-containing protein 3 (RING3-like protein), BRD3, KIAA0043 RING3L
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	BRD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.