



RHEB Antibody, Biotin conjugated

Product Code	CSB-PA624022LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15382
Immunogen	Recombinant Human GTP-binding protein Rheb protein (1-182AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTP-binding protein Rheb (Ras homolog enriched in brain), RHEB, RHEB2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHEB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.