



RIC1 Antibody, HRP conjugated

Product Code	CSB-PA671121LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4ADV7
Immunogen	Recombinant Human RAB6A-GEF complex partner protein 1 protein (951-1233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAB6A-GEF complex partner protein 1 (Connexin-43-interacting protein of 150 kDa) (Protein RIC1 homolog), RIC1, CIP150 KIAA1432
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RIC1