



RIC1 Antibody, Biotin conjugated

Product Code	CSB-PA671121LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q4ADV7
Immunogen	Recombinant Human RAB6A-GEF complex partner protein 1 protein (951-1233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RAB6A-GEF complex partner protein 1 (Connexin-43-interacting protein of 150 kDa) (Protein RIC1 homolog), RIC1, CIP150 KIAA1432
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RIC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.