




Osm Antibody

Product Code	CSB-PA723970LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q65Z15
Immunogen	Recombinant Rat Oncostatin-M protein (26-207AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oncostatin-M (OSM), Osm
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Osm

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA723970LA01RA at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)