



RUNX1T1 Antibody, Biotin conjugated

Product Code	CSB-PA726765HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06455
Immunogen	Recombinant Human Protein CBFA2T1 protein (416-525AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein CBFA2T1 (Cyclin-D-related protein) (Eight twenty one protein) (Protein ETO) (Protein MTG8) (Zinc finger MYND domain-containing protein 2), RUNX1T1, AML1T1 CBFA2T1 CDR ETO MTG8 ZMYND2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RUNX1T1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.