



Thrsp Antibody, HRP conjugated

Product Code	CSB-PA726789HB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q62264
Immunogen	Recombinant Mouse Thyroid hormone-inducible hepatic protein (1-150AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thyroid hormone-inducible hepatic protein (Spot 14 protein) (S14) (SPOT14), Thrsp, S14
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Thrsp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.