



SPOCK1 Antibody, Biotin conjugated

Product Code	CSB-PA730785LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08629
Immunogen	Recombinant Human Testican-1 protein (290-435AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Testican-1 (Protein SPOCK), SPOCK1, SPOCK TIC1 TICN1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SPOCK1