



RABGAP1L Antibody, Biotin conjugated

Product Code	CSB-PA732584LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7ZAP0
Immunogen	Recombinant Human Rab GTPase-activating protein 1-like, isoform 10 protein (11-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab GTPase-activating protein 1-like, isoform 10, RABGAP1L, HHL KIAA0471
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RABGAP1L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.