



RNF20 Antibody, HRP conjugated

Product Code	CSB-PA736198LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VTR2
Immunogen	Recombinant Human E3 ubiquitin-protein ligase BRE1A protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase BRE1A (BRE1-A) (hBRE1) (EC 2.3.2.27) (RING finger protein 20) (RING-type E3 ubiquitin transferase BRE1A), RNF20, BRE1A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RNF20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.