



TENT5C Antibody, Biotin conjugated

Product Code	CSB-PA736203LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VWP2
Immunogen	Recombinant Human Terminal nucleotidyltransferase 5C protein (1-391AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Putative nucleotidyltransferase FAM46C (EC 2.7.7.-), FAM46C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TENT5C