



NCR3LG1 Antibody, Biotin conjugated

Product Code	CSB-PA737873LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q68D85
Immunogen	Recombinant Human Natural cytotoxicity triggering receptor 3 ligand 1 protein (25-262AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Natural cytotoxicity triggering receptor 3 ligand 1 (B7 homolog 6) (B7-H6), NCR3LG1, B7H6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NCR3LG1