




FAT4 Antibody

Product Code	CSB-PA740927LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6V0I7
Immunogen	Recombinant Human Protocadherin Fat 4 protein (4658-4862AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin Fat 4 (hFat4) (Cadherin family member 14) (FAT tumor suppressor homolog 4) (Fat-like cadherin protein FAT-J), FAT4, CDHF14 FATJ
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FAT4

Image



Immunofluorescent analysis of HepG2 cells using CSB-PA740927LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)