



FAT4 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA740927LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q6V0I7 |
| Immunogen | Recombinant Human Protocadherin Fat 4 protein (4658-4862AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Protocadherin Fat 4 (hFat4) (Cadherin family member 14) (FAT tumor suppressor homolog 4) (Fat-like cadherin protein FAT-J), FAT4, CDHF14 FATJ |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | FAT4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |