



MED25 Antibody, Biotin conjugated

Product Code	CSB-PA741144HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q71SY5
Immunogen	Recombinant Human Mediator of RNA polymerase II transcription subunit 25 protein (112-140AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of RNA polymerase II transcription subunit 25 (Activator interaction domain-containing protein 1) (Activator-recruited cofactor 92 kDa component) (ARC92) (Mediator complex subunit 25) (p78), MED25, ACID1 ARC92 PTOV2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MED25