



HIST2H3A Antibody, HRP conjugated

Product Code	CSB-PA744479HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q71DI3
Immunogen	Recombinant Human Histone H3.2 protein (108-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H3.2 (Histone H3/m) (Histone H3/o), HIST2H3A; HIST2H3C; HIST2H3D H3F2 H3FM
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIST2H3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.