




TAF2 Antibody

Product Code	CSB-PA747395LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P1X5
Immunogen	Recombinant Human Transcription initiation factor TFIID subunit 2 protein (400-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription initiation factor TFIID subunit 2 (150 kDa cofactor of initiator function) (RNA polymerase II TBP-associated factor subunit B) (TBP-associated factor 150 kDa) (Transcription initiation factor TFIID 150 kDa subunit) (TAF(II)150) (TAFII-150) (TAFII150), TAF2, CIF150 TAF2B
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TAF2

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA747395LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.